

HERCULES

3000

Pride Mobility Products, Corp.

The Hercules 3000 Lift is the only lift specifically engineered for Pride products and, like the products it is designed to complement, it is strong, durable and very dependable. The Hercules 3000 Lift fits readily into most vehicles and makes transporting your Pride Scooter or Jazzy Power Chair quick and convenient.



Hercules 3000

Designed to lift most of our Pride Jazzy Power Chairs,* Pride Scooters, and PMVs.

FEATURES

- Lifts up to 125 kg (275 lbs.)
- Easily installs in most vehicles
- Conveniently adjustable with quick release pins
- Complete line of docking devices available
- Easily assembles and disassembles
- Our exclusive new pivoting hook for easy maneuvering of Scooter/Jazzy into vehicle
- Transport lock pin added for boom arm stability during transport
- Improved motor gear box design offers smoother & more efficient operation than ever before
- 1-year limited warranty

Various docking devices are available and they are extremely versatile and easily adaptable to virtually any Pride Scooter or Jazzy Power Chair. Attachments are available that will allow you to leave the seat on your Pride Scooter or Jazzy Power Chair depending on the cargo space of your vehicle.*



Pride
Mobility Products
Australia Pty. Ltd.

21 Healey Rd.
Dandenong 3175
Victoria, Australia
phone: (03) 9706-4611
fax: (03) 9706-4622

* Please contact your Pride Provider to confirm Jazzy models/seating systems that the Hercules 3000 is capable of lifting and also to determine if the Hercules 3000 can be fitted in your vehicle.

HERCULES

Pride Mobility Products, Corp.

3000

Technical Information



Straight Boom Arm Dimensions:

Fully Retracted Base and Arm (A)

Total Height: 88 cm (34.75")

Total Length: 77 cm (30.25")

Fully Retracted Base and Fully Extended Arm (B)

Total Height: 99 cm (39")

Total Length: 91 cm (36")

Fully Extended Base and Fully Retracted Arm (C)

Total Height: 102 cm (40")

Total Length: 72 cm (28.5")

Fully Extended Base and Arm (D)

Total Height: 112 cm (44.25")

Total Length: 91 cm (36")

Angled Boom Arm Dimensions:

Fully Retracted Base and Arm (E)

Total Height: 71 cm (28")

Total Length: 76 cm (29.75")

Fully Retracted Base and Fully Extended Arm (F)

Total Height: 95 cm (32.25")

Total Length: 95 cm (37.25")

Fully Extended Base and Fully Retracted Arm (G)

Total Height: 95 cm (37.5")

Total Length: 73 cm (28.75")

Fully Extended Base and Arm (H)

Total Height: 95 cm (37.5")

Total Length: 95 cm (37.25")