

**SPECIFICATIONS**

	<b>C1</b>	<b>C5</b>	<b>C6</b>	<b>C-101</b>	<b>LC-101 Euro Leather</b>	<b>LC-900 Media Chair</b>
<b>Back Width:</b>	51 cm (20")	71 cm (28")	71 cm (28")	72 cm (28.5")	58 cm (23")	79 cm (31")
<b>Seat-To-Floor:</b>	41 - 46 cm (16" - 18")	46 - 51 cm (18 - 20")	48 cm (19")	47 cm (18.5")	47 cm (18.5")	52 cm (20.5")
<b>Seat Depth:</b>	51 cm (20")	46 - 51 cm (18 - 20")	53 cm (21")	51 cm (20")	51 cm (20")	56 cm (22")
<b>Seat Width:</b>	46 cm (18")	56 cm (22")	56 cm (22")	51 cm (20")	51 cm (20") <sup>1</sup>	56 cm (22")
<b>Top of Back to Seat:</b>	66 cm (26")	66 - 71 cm (26 - 28")	66 - 71 cm (26-28")	70 cm (27.5")	70 cm (27.5")	58.5 cm (23")
<b>Weight Capacity:</b>	113 kg (250 lbs.)	170 kg (375 lbs.)	136 kg (300 lbs.)	147 kg (325 lbs.)	147 kg (325 lbs.)	170 kg (375 lbs.)
<b>Colour Options:</b>	Arctic Blue Chocolate Oatmeal	Arctic Blue Chocolate Oatmeal	Arctic Blue Chocolate Oatmeal Royal Blue	Cocoa Fern Sandal	<b>Leather:</b> Beige Brown Black Red	Black (Lexis Vinyl)

  

	<b>D30</b>	<b>560</b>	<b>660</b>	<b>LL-805</b>	<b>T-3</b>	<b>LC-440 Power Stand Assist</b>
<b>Back Width:</b>	74 cm (29")	71 cm (28")	71 cm (28")	77.5 cm (30.5")	69 cm (27")	53.5 cm (21")
<b>Seat-To-Floor:</b>	48 cm (19")	48 cm (19")	48 cm (19")	48 cm (19")	48 cm (19")	49 cm (19")
<b>Seat Depth:</b>	56 cm (22")	48 cm (19")	51 cm (20")	51 cm (20")	51 cm (20")	46 cm (18")
<b>Seat Width:</b>	56 cm (22")	48 cm (19")	51 cm (20")	51 cm (20")	51 cm (20")	54.5 cm (21.5")
<b>Top of Back to Seat:</b>	66 cm (26")	63.5 cm (25")	63.5 cm (25")	68.5 cm (27")	68.5 cm (27")	69 cm (27")
<b>Weight Capacity:</b>	170 kg (375 lbs.)	170 kg (375 lbs.)	170 kg (375 lbs.)	136 kg (300 lbs.)	170 kg (375 lbs.)	136 kg (300 lbs.)
<b>Colour Options:</b>	Ocean Burgundy Spruce Latte <b>Levante Vinyl:</b> Burgundy Royal Blue Pine Mountain	Ocean Burgundy Spruce	Latte	Cocoa Fern Sandal	Black (Leather) Oatmeal	Cocoa Fern Sandal

**BENEFITS**

- Engineered furniture laminate/hardwood frame (C5 & C6 both have steel frames)
- Patented, quiet & smooth lift system
- Large side pocket\*
- Hand control can be easily moved to left-hand side
- Integrated emergency battery backup (batteries not included)
- External transformer for easy access to battery backup
- 2-year limited warranty

Note: Measurements may vary +/- 2.5 cm (1"). Specifications subject to change without prior notice.

\*Large side pocket not available on LC-900 or LC-440.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



**Pride® Lift Chairs**  
Experience the Luxury

# Pride Lift Chairs

*Experience the Luxury*

Pride® Lift Chairs are designed to be the ultimate in lift chair comfort, style and performance. Available in several beautiful fabrics and styles to fit any décor, Pride Lift Chairs also feature a 2-year limited warranty. The state-of-the-art lift system utilized in Pride Lift Chairs will recline you further back than any other lift chair. Plus, Pride Lift Chairs are equipped with an emergency battery backup system should the power fail in your home.

## LIFT CHAIR DESCRIPTIONS

**C1**  
3-position, petite to small lift chair with fully removable upholstery components

**C5**  
3-position, medium to large chaise lounger with fully removable upholstery components

**C6**  
Infinite-position, medium to large chaise lounger with fully removable upholstery components and dual motors which enable independent recline/footrest function

**C-101**  
3-position, medium to large lift chair

**LC-101 Euro Leather**  
3-position, medium to large lift chair

**LC-900 Media Chair**  
3-position, medium to large chaise lounger, available with optional media wedge<sup>2</sup>

**D30**  
3-position, medium to large chaise lounger

**560**  
3-position, small to medium chaise lounger

**660**  
Infinite-position, small to medium chaise lounger with dual motors that enable independent recline/footrest function

**LL-805**  
2-position, medium to large wall hugger chaise lounger which requires only 10 cm (4") of space between wall and chair

**T-3**  
Infinite-position, medium to large lift chair featuring Trendelenburg position (elevates legs above head for greater relaxation) and zero gravity position

**LC-440 Power Stand Assist<sup>3</sup>**  
Power seat lift

## AVAILABLE FABRICS <sup>1</sup>



<sup>1</sup> Not all fabrics are available for every model. Please see Specifications for details.  
<sup>2</sup> Media wedge sold separately and not included with purchase of the LC-900.  
<sup>3</sup> This chair does not recline.

## LIFT CHAIR MODELS



**C1**

Accented Pillow-back Design  
Shown in Arctic Blue



**C5**

Sewn Pillow-back Design  
Shown in Arctic Blue



**C6**

Sewn Pillow-back Design  
Shown in Chocolate



**C-101**

Button-back Design  
Shown in Sandal



**LC-101 Euro Leather**

Split-back Design  
Shown in Beige Leather



**LC-900 Media Chair<sup>2</sup>**

Sewn Pillow-Back Design  
Shown in Black (Lexis Vinyl)



**D30**

Split-back Design  
Shown in Burgundy



**560**

Pillow-back Design  
Shown in Spruce



**660**

Pillow-back Design  
Shown in Latte



**805**

Sewn Split-back Design  
Shown in Cocoa



**T-3**

Mongolian-back Design  
Shown in Oatmeal



**LC-440 Power Stand Assist<sup>3</sup>**

Shown in Sandal



**PLEASE NOTE:** 2-position is a TV recline.  
 3-position is a full recline.  
 Infinite-position is a full recline with independent leg rest function.